

Documentation for EMSAR

version 2.0.0

Author : Soo Lee (duplexa@gmail.com)

Author Affiliation : Center for Biomedical Informatics, Harvard Medical School, Boston, MA 02135, USA

Overview

EMSAR is a C program that quantifies gene/transcript expression from RNA-seq data. It first builds an rsh index from a transcriptome fasta file given a specified library type and read length (eg. Paired-end, unstranded, 101 bp), then uses it for speedy processing of alignments to estimate the expression levels. The rsh index contains information of sequence sharing among different genes and mRNA isoforms.

EMSAR consists of two parts : emsar-build and emsar. The first builds an rsh index. The second takes alignment files and produce expression estimates. Alternatively, the user can also build an rsh index internally by starting with a transcriptome fasta file when running emsar on RNA-seq read alignments. In the latter case, read length is automatically obtained from the first bam file and the user can choose either to create an rsh index file or to use the information only internally. It is highly recommended to use emsar-build to create a reusable rsh index file, which greatly saves time at later steps. It will not affect the actual quantification result (unless the user specifies incorrect read length or library type when building the rsh index).

Building an rsh index (emsar-build)

To build an rsh file, the user must specify a transcriptome fasta file. Note that this fasta file must be identical to the one used for mapping RNA-seq reads. Otherwise, the final estimation accuracy may be lower. EMSAR will not check the identity of the two transcriptomes and it is the user's responsibility to ensure this.

The read length (range) must also be specified. For a paired-end library, the read length must be identical for all RNA-seq data and it should match the one specified for the rsh index. For a single-end library, a range of read lengths can be set.

It may take several hours for a complex transcriptome such as a recent human ENSEMBL annotation, particularly for a paired-end library. It is recommended to use multi-threading. The run time is proportional to the range of fragment lengths to cover. By default, it is 1-400 bp, but if the user wants to ensure that it can be applied to RNA-seq data with longer fragment lengths, there is an option to change this range, at the cost of speed. RNA-seq fragments beyond the specified range will not be included for calculation at a later step.

Additionally, the user must specify an output directory and an output file prefix.

The following shows an example run.

```
emsar-build -P myTranscriptome.fa 101 myOutdir myOutprefix
```

This command will create an rsh file, `myOutdir/myOutprefix.rsh`, for paired-end (-P), unstranded (default) 101 bp x 2 reads, covering fragment lengths of 1-400 bp (default). An example `myOutprefix` could be `myTranscriptome.L101x2`.

Multi-threading option

It is recommended that the user take advantage of the multi-threading option for speedy computation. Multi-threading can be turned on with a -p option. An example run using 8 threads would look like:

```
emsar-build -P -B -p 8 myTranscriptome.fa 101 myOutdir myOutprefix
```

Main program (emsar)

EMSAR requires an rsh index file and an alignment file (bam, sam or bowtie1 default output format) as input. It assumes that the alignment file was on the same transcriptome fasta file (not on the genome). If the two do not match, the program may produce an unpredictable outcome. Likewise, the read length and library type (paired-end vs single-end, unstranded vs forward-stranded vs reverse-stranded) must match. Alternatively, the transcriptome fasta file can be fed directly (-x) instead of an rsh index file (-I), but it will take much longer.

Additionally, the user must specify an output directory and an output file prefix.

A simple example run would look as follows:

```
emsar -P -B -I myTranscriptome.rsh myOutdir myOutprefix myPESample.bam  
emsar -P -B -x myTranscriptome.fa myOutdir myOutprefix myPESample.bam
```

Here, myOutdir and myOutprefix must represent the corresponding samples.

Note that -P option is used for a paired-end sample and -B option is used for a bam file. By default, single-end, unstranded library is assumed and bowtie1 output format is a default alignment format. For a sam file with a single-end library, EMSAR could be run as follows:

```
emsar -S -I myTranscriptome.rsh myOutdir myOutprefix mySESample.SAM
```

Multi-threading option

It is recommended that the user take advantage of the multi-threading option for speedy computation. Multi-threading can be turned on with a -p option. An example run using 8 threads would look like:

```
emsar -P -B -p 8 -I myTranscriptome.L75x2.rsh myOutdir myOutprefix  
myPESample.bam
```

Alignment

Theoretically, any transcriptome-aligned alignment file can be used, but for a good performance, it is recommended to use the following parameters with bowtie1.

```
-a -v 2 -X 1000 -m 100
```

The -a and -v 2 options are recommended since EMSAR internally filters all best alignments, assuming that all alignments up to a certain number of mismatches are reported. For this filtering, EMSAR assumes an auxiliary field in the bam or sam file (the MD field) that reports the number of mismatches. Bowtie1 default output format contains this information, which is again used by EMSAR.

The -X option specifies the maximum fragment length in bp for paired-end data. It is usually safer to set this number to as high as 1000 even if the actual fragment size is smaller than that. Note that EMSAR will use the fragments within the range used to build the rsh index. If the fragment lengths tend to be larger, then one can use the -F (--maxfragen) option of either emsar-build or emsar to increase the maximum fragment size. For single-end data, fragment length-related parameters can be ignored.

The -m 100 option suppresses alignments if more than 100 alignments exist for a read (pair). EMSAR internally filters out alignments similarly based on the -k option of emsar-build and emsar, which is set to be 100 by default. It is not necessary to use the -m option for bowtie run, but if one wants to use it, it would be safe to set this number to be at least equal to or larger than the -k value set by emsar-build and emsar.

Note on the alignment with mismatches

There is a chance that reads aligned with mismatches are not well matched with the computed length assuming no mismatch. This problem is relatively well handled with EMSAR's internal alignment filtering step, but under rare cases, it is possible that a read is mapped to two locations, each with one mismatch on a different nucleotide. These two locations will not have the identical sequences, therefore they will not be considered as a proper segment with a nonzero length. These rare alignments that are in conflict with the segments defined on a perfect-match scenario are excluded from the calculation.

Transcriptome fasta file

It is strongly recommended not to use a polyA-attached transcriptome fasta file. This causes very slow suffix array sorting and sometimes cause an unknown (therefore unfixed) problem. The fasta file may contain duplicated sequences, but in that case, the user must be cautious in interpreting the final fpkm values.

Streaming of alignments

An alignment file can be piped directly from an alignment program. For example, if the user begins with fastq files and aligned them to a transcriptome index using bowtie1, the command would look like:

```
bowtie -a -v 2 -m 100 myTranscriptomeIndex mySample.fastq | emsar -I  
myTranscriptome.rsh myOutdir myOutprefix
```

Note that the transcriptome index file must have been built from the same transcriptome fasta file. The above command assumes a single-end library. To stream as a sam or bam file, one could use the following command:

SAM :

```
bowtie -a -v 2 -m 100 -S myTranscriptomeIndex mySample.fastq | emsar  
-S -I myTranscriptome.rsh myOutdir myOutprefix
```

BAM :

```
bowtie -a -v 2 -m 100 -S myTranscriptomeIndex mySample.fastq |  
samtools view -hbS - | emsar -B -I myTranscriptome.rsh myOutdir  
myOutprefix
```

Multi-sample option

A user can specify a text file containing a list of alignment files instead of a single alignment file, along with the -M option (--multisample). For example, the following command would run EMSAR in a multi-sample mode.

```
emsar -P -B -M -I myTranscriptome.rsh myOutdir myOutprefix
```

```
myBam.list.txt
```

myBam.list.txt could look like:

```
myInputdir/myPESample1.bam
myInputdir/myPESample2.bam
myInputdir2/myPESample3.bam
```

Note that the -P and -B option is used to match the alignment file type in the alignment list file. Currently, a multi-sample run with mixed library types or alignment file formats is not supported. Alignment streaming is not supported in the multi-sample mode.

Library type

The strandedness of the library can be specified with the -s (--strand_type) option. By default, 'ns' (unstranded) is used for both single-end or paired-end libraries. For a stranded single-end library, one can use 'ssf' (forward stranded) or 'ssr' (reverse stranded). The former assumes that the reads are mapped to the sense strands of the transcriptome. The latter assumes antisense. For a paired-end library, 'ssfr' or 'ssrf' can be used. The former assumes that the first mate is forward and the second mate is reverse. 'ssrf' is the opposite. An example usage is provided below:

```
emsar -P -B -s ssfr -I myTranscriptome.rsh myOutdir myOutprefix
myStrandedPESample.bam
```

Maximum number of alignments per read

The -k (--max_repeat) option of EMSAR sets the maximum number of alignments per read. Any read with more than this number will be ignored for both sample read counts and for computing virtual (effective) lengths. If one uses bowtie for alignment, the -m option can be used to restrict the number of alignments per read. If this is set to be smaller than the -k option of EMSAR, the accuracy may be slightly sacrificed because the maximum number of alignments for computing virtual length is different from that of the sample. Therefore, it is recommended that the user either matches the -m option of bowtie with -k option of EMSAR or make sure that the -m option is set to be larger than the -k option. A similar recommendation would also apply for other aligners.

If, for example, the user wants to use reads mapped to 3000 times, then, the bowtie and EMSAR commands could look like:

```
bowtie -S -a -v 2 -m 3000 -X 1000 myTranscriptomeIndex -1
mySample.R1.fastq -2 mySample.R2.fastq | samtools view -hbS - >
mySample.bam
emsar-build -P -k 3000 myTranscriptome.fa myOutdir myOutprefix
emsar -P -B -k 3000 -I myTranscriptome.rsh myOutdir myOutprefix
mySample.bam
```

Segment information

The list of all segments with their associated sequence-sharing set and transcripts(separated by '+') is generated with the -g option. The output .segments file also contains the effective length of each segment for each transcript (separated by comma), observed read count and expected read count for each segment. The expected read count is based on the estimated FPKM values and the effective lengths. An example .segments file looks like :

segment_id	sequence	sharing_set_id	transcript_ids	transcript_names	eff.length	Readcount	expected_Readcount
c0	s0	t0	ENSVPAT00000007931	474.005893 39 39.000000			
c1	s1	t1	ENSVPAT00000004546	0.000000 0 0.000000			
c2	s2	t2	ENSVPAT00000003367	1720.000000 84 84.000000			
c3	s3	t3	ENSVPAT00000009296	608.011785 59 59.000000			
....							
c12916	s2644	t4871,t7233,t8590	ENSVPAT00000004400+ENSVPAT00000000132+ENSVPAT00000010274	20.000000 0 4.250379			
c12917	s2644	t4871,t8238,t8590	ENSVPAT00000004400+ENSVPAT00000004654+ENSVPAT00000010274	1.639434 1 0.306191			
c12918	s307	t4892,t7897,t11970	ENSVPAT000000004120+ENSVPAT00000008179+ENSVPAT00000001600	5.005893 0 0.950490			

A segment with a single transcript (eg. the first line, segment_id c0) represents the region unique to the transcript. The effective length (eff.length) is the length of the unique region of this transcript, as the number of possible alignments adjusted for fragment length distribution. A segment associated with multiple transcripts (eg. The last line, segment_id c12918) represents the region shared by (and only by) the listed transcripts and the eff.length is the length of such shared region. Readcount and expected_Readcount represent the observed and expected read counts for such unique or shared region.

An example command would look like:

```
emsar -P -B -g -I myTranscriptome.rsh myOutdir myOutprefix
mySample.bam
```

Output logging

A updating progress bar is printed to standard output when the -v (verbose) option is used. It is highly recommended not to use this option when standard output is logged. By default, major steps are logged with time stamps. Logging can be turned off by using the -q option.

utils

The following utils are provided as a part of the EMSAR package.

FPKM2gFPKM.pl
transcript_stats.3.pl
post_processing.pl

FPKM2gFPKM.pl

This script sums the isoform-level estimates to obtain gene-level expression level estimates. Any gene-to-isoform mapping can be used. For example, one can use a tab-delimited text file containing an ENSEMBL gene ID and an ENSEMBL transcript ID per line (.g2t file) or a file containing a gene symbol and an ENSEMBL transcript ID per line (.s2t file). The transcript IDs are the same as the ones used in the fasta file and the gene IDs (symbols) are what the user wants the expression levels to be computed for. The example files are shown below.

An example ENSEMBL .g2t file:

```
ENSG00000213366 ENST00000369829
ENSG00000224246 ENST00000448586
ENSG00000034677 ENST00000522182
```

```

ENSG00000160323  ENST00000371910
ENSG00000136883  ENST00000491059
ENSG00000054690  ENST00000559981
ENSG00000258905  ENST00000553993

```

An example ENSEMBL .s2t file:

```

GSTM2  ENST00000369829
SFTA2  ENST00000448586
RNF19A ENST00000522182
ADAMTS13      ENST00000371910
KIF12  ENST00000491059

```

An example usage would look like:

```
FPKM2gFPKM.pl myTranscript.g2t mySample.fpkm > mySample.gfpkm
```

The output file of the script would look as below:

geneID	FPKM	iReadcount	iReadcount.int	TPM
RP11-560A15.3	0	0	0	0
RPS11	1168.165255	19756	19756	2114.136652
CREB3L1	90.215036	9361.145276	9361	163.270489
RPL10P14	0	0	0	0

transcript_stats.3.pl

usage : transcript_stats.3.pl fastafile segmentfile g2tfile

The script generates a series of traits associated with each transcript, including GC content, number of isoforms, effective length unique to a transcript ('isoform_unique_length'), effective length unique to a gene ('gene_unique_length'), effective length shared with other genes ('multi_gene_length') and their proportions. In order to use this script, one should have run EMSAR with the -g option, to first create a .segments file which is used as one of the input files here. The other two input files are the transcriptome fasta file and the g2t file as described above. An example output file is shown below (with line numbers to help tracking lines):

1:	transcript_id	gene	transcript_length	GC_content	nIsoforms		
	total_effective_length	isoform_unique_length	gene_unique_length	multi_gene_length			
	gene_unique_isoform_common_length	isoform_unique_proportion	gene_unique_proportion				
	gene_unique_isoform_common_proportion						
2:	ENST00000369829	GSTM2	914	0.509	15	753.258783	322.000000 752.626423
	0.63236	0	0.427	0.999	0		
3:	ENST00000598454	LAMTOR5-AS1	711	0.495	13	551.023467	373.157403
	549.943391	1.080076	0	0.677	0.998	0	
4:	ENST00000448586	SFTA2	451	0.556	12	0	0 0 0 NA
	NA	NA					
5:	ENST00000522182	RNF19A	544	0.522	15	384.023466	165.269720 384.023466
	0	0	0.430	1	0		
6:	ENST00000371910	ADAMTS13	1074	0.583	17	914.023467	151.000000
	914.023467	0	0	0.165	1	0	

An example usage would look like:

```
transcript_stats.3.pl myTranscriptome.fa mySample.segments
myTranscript.g2t > mySample.stats
```

post_processing.pl

usage : post_processing.pl fpkm_file_dir g2t_file

The script converts all the isoform-level fpkm files (direct output of EMSAR) in a specified directory to gene-level fpkm (gfpkm) files, given a g2t file described above, merges the gene-level read counts and TPM values of multiple samples into two big tables (gReadcount.all and gTPM.all).